

B007 : Asphalt Splitter

This instrument is used to break and crumble asphalt samples to facilitate bulk density tests and laboratory testing purposes by reducing the granulating time in few minutes with high quality results.

A rough asphalt sample of approx. 1 kg is poured into the stainless steel bowl equipped, at its bottom, with three rotating paddles. The cover is closed, the machine started, and the three paddles break and crumble the sample in grains.

The bowl is now tilted to discharge the sample into the self supporting suitable pan supplied with.

A hinged cabinet reducing noise protects the bowl. When opening the cover while the splitter is working, a microswitch automatically stops the machine. Conforming to CE Safety Directive.

Power supply: 230V 50Hz 1500W

Dimensions: 490x730x855 mm

Weight: 85 kg approx.

Spare part :

B007-11 Kit of paddles



B007

BAG-MacBen BV
Oostjachtpark 12,
B-9100 Sint-Niklaas

T: +32 3 771 48 04 info@macben.eu
www.macben.eu

BAG-MacBen
Smederijstraat 2,
NL-4814 DB Breda

T: +31 765 30 23 73 sales@macben.eu
www.macben.eu